

BCAL gain Calibration using π^0 s



Karthik Suresh



University
of Regina



Faculty of
Science





NSERC
CRSNG

Outline





1. Detectors in GlueX
2. BCAL- design and working principle
3. BCAL Gain Calibrations
 - a. Energy determination algorithm
 - b. Calibration using π^0
 - c. To address the presence of nonlinearity after the calibration.



- BCAL detects $\rightarrow n, \gamma, p, \pi^+, \pi^-$
- BCAL along with tracking chambers gives the PID information

Barrel Calorimeter (BCAL)




- ❑ 4 m long detector with $r_{<} = 65 \text{ cm}$ and $r_{>} = 90 \text{ cm}$
- ❑ Polar acceptance $12^\circ < \theta < 126^\circ$
- ❑ Segmented into 48 azimuthal modules with each module containing 40 SiPMs. Provides almost 2π coverage
- ❑ 1 module has 185 layers of Pb sheet with 184 layers of double clad scintillating fibres.
- ❑ Total 685000 fibres connected to 3840 SiPMs (Resistant up to 2T) with 1536 readouts


BCAL (continued)



BCAL before inserting into the bore magnet



SiPM attached readout module




Summing scheme


$$\frac{\sigma(\mathbf{E})}{\sqrt{E}} = \frac{5.2\%}{\sqrt{E}} \oplus 3.6\%$$

$$\sigma(\mathbf{z}) = 3\text{cm} \quad \sigma(t) \approx 200\text{ps} @ 1\text{GeV}$$

Determination of energy and calibration using π^0 in BCAL




- Energy of photons are determined by the SiPM's attached at the end of the Sci-Fibres.
- The position of photon shower is reconstructed using the timing information. The vertex of reaction is calculated from charged particle.
- π^0 s offer large statistics in the energy range of interest (0.5-3 GeV)





The Calibration Algorithm

1. Events with $\pi^0 \rightarrow 2\gamma$ with at least 2 charged particles.
2. Determine photon energies ($E_{\gamma_1}, E_{\gamma_2}$)
3. Find event vertex using the charged particles.
4. Select channels in layer 1 and 2 having more than 50% of the deposited energy (E_{γ_1} and E_{γ_2})
5. Reconstruct invariant mass of π^0 using the E_{γ} 's and the vertex^[1], and fit it. $m^2 = 2E_1E_2(1 - \cos(\psi))$
6. Form ratio with PDG invariant mass (0.135 GeV) to the measured mean of invariant mass and use to adjust gains for each channel
7. The entire procedure is iterated till tolerance is less than 1%






[1] <https://doi.org/10.1016/j.nima.2018.04.006>

Calibration plots




Calibration plots (continued)




Checks in Different energy ranges


σ/μ vs. Iteration (E1, E2 > 500 MeV)




σ/μ vs. Iteration (E1, E2 > 700 MeV)



σ/μ vs. Iteration (E1, E2 > 900 MeV)




σ/μ vs. Iteration (E1, E2 > 1100 MeV)




Nonlinearity after the calibration

- SiPM and FADC Saturation: FADC \rightarrow Constant pulse shape normalisation: SiPM No correction
- An empirical nonlinear function is fit to the plot and is used to correct for nonlinearity

π^0 Mean vs Energy bin



π^0 Mean vs Energy bin



Conclusions

- ❑ BCAL is a key detector in the GlueX experiment.
- ❑ The BCAL performance is monitored online.
- ❑ A new gain calibration for the BCAL is needed for every run period BCAL is necessary. The π^0 method shown here works well.
- ❑ SiPM saturation at high energies have also been employed.

References

- T. Beattie et al , NIMA 896 (2018) 24-42 <https://doi.org/10.1016/j.nima.2018.04.006>
- E.Smith “Note on Saturation in BCAL SiPMs” [GlueX-doc-3737-v4](#)

Thank you



GlueX Collaboration Meeting

May 15-17, 2017

Jefferson Lab • Newport News, VA